
# Poseidon






## Hydrophilic Nitinol Guide Wire

EN	<i>English</i>	<b>Hydrophilic Nitinol Guide Wire</b>
BG	<i>Bulgarian/Български език</i>	<b>Хидрофилен нитинолов водач</b>
HR	<i>Croatian/Hrvatski</i>	<b>Hidrofilna žica vodilica od nitanola</b>
CS	<i>Czech/Česky</i>	<b>Hydrofilní nitanolový vodící drát</b>
DA	<i>Danish/Dansk</i>	<b>Hydrofil nitanolguidewire</b>
NL	<i>Dutch/Nederlands</i>	<b>Hydrofele Nitinol voerdraad</b>
ET	<i>Estonian/Eesti</i>	<b>Hüdrofilne nitanoolist juhtetraat</b>
FI	<i>Finnish/Suomi</i>	<b>Hydrofilinen nitanolihjoinlanka</b>
FR	<i>French/Français</i>	<b>Fil-guide Nitinol hydrophile</b>
DE	<i>German/Deutsch</i>	<b>Hydrophiler Nitinol Führungsdraht</b>
EL	<i>Greek/Ελληνικά</i>	<b>Υδροφιλο οδηγό σύρμα ντινιόλης</b>
HU	<i>Hungarian/Magyar</i>	<b>Hidrofil nitanol vezetődrót</b>
IT	<i>Italian/Italiano</i>	<b>Filo guida in nitanol idrofilo</b>
LV	<i>Latvian/Latviešu</i>	<b>Nitanola vadītājstīga ar hidrofilu pārklājumu</b>
LT	<i>Lithuanian/Lietuvių k.</i>	<b>Hidrofilinė nitanolo kreipiamoji viela</b>
NO	<i>Norwegian/Norsk</i>	<b>Hydrofil nitanolledsonde</b>
PL	<i>Polish/Polski</i>	<b>Hydrofilowy przewodnik nitanolowy</b>
PT	<i>Portuguese/Potuguês</i>	<b>Fio-guia de nitanol hidrofílico</b>
RO	<i>Romanian/Română</i>	<b>Fir de ghidaj hidrofil din nitanol</b>
SK	<i>Slovak/Slovenčina</i>	<b>Hydrofil nitanol vodiaci drôt</b>
SL	<i>Slovenian/Slovenščina</i>	<b>Hidrofilna nitanolna vodilna žica</b>
ES	<i>Spanish/Español</i>	<b>Alambre guía hidrófilo de nitanol</b>
SV	<i>Swedish/Svenska</i>	<b>Hydrofil nitanolledare</b>
TR	<i>Turkish/Türkçe</i>	<b>Hidrofilik Nitinol Kılavuz Tel</b>





**EN: Explanation of symbols used on the package labels - BG: Значение на символите върху етикетите на опаковката - HR: Objašnjenja simbola koji se koriste na naljepnici na pakiranju: - CS: Vysvětlení symbolů použitých na štítcích obalů - DA: Forklaring af symboler anvendt på emballagens etiketter**

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## Poseidon Hydrophilic Nitinol Guide Wire

### Description:

The Poseidon Hydrophilic Nitinol Guide Wire is EO sterilized, non-pyrogenic and available in the following variants:

Shaft stiffness	Tip stiffness	Tip subtype (shape)	Available variants
<ul style="list-style-type: none"> <li>Standard</li> <li>Stiff</li> </ul>	<ul style="list-style-type: none"> <li>Floppy</li> <li>Soft</li> <li>Standard</li> </ul>	<ul style="list-style-type: none"> <li>Straight</li> <li>Angled</li> <li>Double ended</li> <li>Reshapable</li> </ul>	<ul style="list-style-type: none"> <li>Guide wire diameter</li> <li>Tip shape</li> <li>Tip stiffness</li> <li>Shaft stiffness</li> <li>Guide wire length</li> </ul>

For further information, see SP Medical's website.

### Intended purpose:

The Poseidon Hydrophilic Nitinol Guide Wire is indicated for general intravascular and coronary arterial use to aid in the selective placement of interventional devices during diagnostic and / or therapeutic procedures.

### Contraindication:

**Not for use in the cerebral vasculature.**

### Warnings:

1. SINGLE USE. This guide wire is intended for single use only. Reusing the guide wire involves high risk of contamination and locking of the guide wire inside the interventional device due to wear-out of the hydrophilic coating.
2. This guide wire is provided in sterile condition. Do not re-sterilize or re-use the guide wire.
3. Do not withdraw through a metal needle.
4. Avoid interaction with metal / hard plastic devices with sharp edges or burrs that may separate the coating or polymer from the guide wire.
5. Never advance or withdraw the guide wire against resistance until the cause of the resistance is determined by fluoroscopy.
6. Do not attempt to move the guide wire without observing the resulting tip response.
7. Do not use alcohol, antiseptic solutions, other solvents or dry gauze on the guide wire, as this may damage the hydrophilic coating.
8. The guide wire contains a metallic core. Do not use with any inappropriate equipment (e.g., MRI).
9. To be used before the expiry date stated on the package.

### Precautions:

1. Prior to opening, the sterile package should be carefully inspected to see if it is still intact. The inspection must be carried out just before the guide wire is used and in good light conditions. Sealing quality on the transparent side of the peel-pouch must be carefully inspected for channels and cracks. Both front and back side must be inspected for pinholes or other packaging defects such as abrasions or cracks. Use a magnification tool if in doubt of sterile barrier integrity.
2. Do not use if the package / sterile barrier is broken / damaged.
3. Prior to use, carefully inspect the guide wire for damages.
4. Do not use damaged guide wires.
5. SP Medical does not recommend a particular technique or device configuration for the use, and the steps contained in these instructions for use are for information purposes only.
6. The guide wire should only be selected and used by experienced physicians, trained in invasive techniques, the use of guide wires and familiar with side effects and hazards commonly associated with interventional procedures.
7. The physician should evaluate the appropriateness according to individual patient condition and his or her medical training and experience.
8. It is always the physician's responsibility to determine and ensure the patient's suitability for the procedure where the guide wire is used.

### Potential side effects:

Possible complications include but are not limited to the following:

- Vessel / tract trauma
- Vessel / tract narrowing / occlusion
- Infection and inflammation
- Allergic reaction

### Serious incident:

Every serious incident that has occurred in relation to the device should be reported without delay to the manufacturer and the competent authority of the Member State in which the user and / or patient is established. For further information, please contact the manufacturer.

### Compatibility:

The guide wire may be used in combination with interventional devices in a CathLab environment. The end-user is solely responsible for selection of the appropriate device. There should be at least 0.0004" (0.01 mm) clearance between the lumen of the interventional device and the guide wire.

### Preparations before use:

1. Before removing the guide wire from the dispenser, use a 20 ml syringe to fill the dispenser with saline (sterile) and allow the guide wire to hydrate for at least 30 seconds.
2. Carefully remove the guide wire from the dispenser.
3. Inspect the guide wire thoroughly to make sure that it is not damaged.
4. Confirm the compatibility of the guide wire diameter with the interventional device before use.

### Directions for use:

1. When introducing the guide wire into the interventional device, ensure that at least 2 centimeters of the guide wire extend from the proximal hub. This will prevent the guide wire from slipping inside e.g., the catheter.
2. To aid the selective placement of the interventional device into a particular vessel, gently rotate the proximal end of the guide wire as it is advanced forward.
3. To prevent contrast agent crystallization / clotting, a continuous saline flush should be maintained between the interventional device and the guide wire during the procedure. The size of the syringe used to flush the lumen should be adapted to the length and diameter of the interventional device.
4. Between uses during the same procedure, place the guide wire in a saline-filled container, or fill the dispenser with saline and reinsert the guide wire in the dispenser, distal end first. Avoid wiping with damp cloths as particles from the cloth can adhere to the surface of the guide wire. Make sure to leave a part of the proximal end outside the dispenser.
5. When the procedure is finished, remove the guide wire and interventional device according to procedural protocol.
6. Discard the product after use according to local instructions for hazardous waste.

### Storage:

To be stored under cool, dark and dry conditions.

Humidity limits: 30% - 70%. Temperature limits: 2°C - 49°C

### Expiry date:

The guide wire must be used no later than stated on the label.

### Disposal:

Disposal shall be handled in accordance with accepted medical practice and applicable local laws and regulations for hazardous waste.

### Safety information:

The IFU is delivered in hard copy together with the products and is kept updated and available for electronic download on SP Medical's website.

Summary of safety and clinical performance can be obtained via Eudamed by entering SP Medical as legal manufacturer and the basic UDI-DI: 5710086-GWPOLZ.

<https://ec.europa.eu/tools/eudamed/#/screen/search-device>

### Disclaimer of liability:

SP Medical A/S is not liable for defects / deterioration resulting from abnormal use or modifications made to the product. Under these circumstances, the product is not covered by the guarantee. SP Medical A/S disclaims liability for direct or indirect injuries that may occur as a consequence of the product being modified or wrongly used.

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