

**Movable Guide Wire**



**Dimensions**

Ø 0.032" / 0,81 mm  
Ø 0.035" / 0,89 mm  
Ø 0.038" / 0,96 mm

**Length**

100 – 450 cm

**Surface**

Stainless steel  
PTFE coated

**Distal tip form**

Straight tip  
J-tip radius: 1,5 – 3 – 6 – 15 mm


**Shaft stiffness**

Standard  
Heavy duty  
Customised

**Mode of delivery**

Sterile – CE marked  
Bulk OEM

**Customised products on request**



SP Medical develops, manufactures and provides plastics and surface solutions to the global medical devices and pharmaceutical industries. Our core capabilities are injection molding, surface treatment and processing of plastic products. At our state-of-the-art facilities in Denmark and Poland, we manufacture components, medical devices, single-use devices, and tablets with active pharmaceutical ingredients, and we apply function-enhancing coatings to a number of products.

SP Medical is a company in SP Group A/S, an international plastics group headquartered in Denmark. SP Group is listed on NASDAQ OMX Copenhagen.

[www.sp-medical.com](http://www.sp-medical.com)