


CHAMPION PTCA Guide Wire – Heavy duty tip for highly stenosed lesions



Dimension

Ø 0.014" / 0,36 mm

Length

195 cm

300 cm

Distal tip form

Straight

J-tip

Stent zone stiffness

Standard

Tip stiffness (Gram Force)

3,0gf

4,5gf

6,0gf

9,0gf

12,0gf

Radiopaque length

3 cm

Coil surface

Uncoated – Stainless Steel

Hydrophilic coated

Shaft surface

PTFE coated

Mode of delivery

Sterile – CE marked

Bulk OEM



SP Medical develops, manufactures and provides plastics and surface solutions to the global medical devices and pharmaceutical industries. Our core capabilities are injection molding, surface treatment and processing of plastic products. At our state-of-the-art facilities in Denmark and Poland, we manufacture components, medical devices, single-use devices, and tablets with active pharmaceutical ingredients, and we apply function-enhancing coatings to a number of products.

SP Medical is a company in SP Group A/S, an international plastics group headquartered in Denmark. SP Group is listed on NASDAQ OMX Copenhagen.

www.sp-medical.com