

Schüller Guide Wire



Dimensions

Ø 0.025" / 0,63 mm
Ø 0.035" / 0,89 mm

Length

70 cm
80 cm
90 cm
100 cm
110 cm
150 cm

Surface

Stainless steel
PTFE coated

Distal tip form

Straight tip
J-tip radius: 1,5 – 3 – 6 – 15 mm

Distal flexible length

5 cm
7 cm
10 cm


Proximal flexible length

30 cm
33 cm
35 cm
38 cm
40 cm
43 cm
45 cm
50 cm
53 cm
55 cm
60 cm
63 cm
65 cm

Mode of delivery

Sterile – CE marked
Bulk OEM

Customised products on request



SP Medical develops, manufactures and provides plastics and surface solutions to the global medical devices and pharmaceutical industries. Our core capabilities are injection molding, surface treatment and processing of plastic products. At our state-of-the-art facilities in Denmark and Poland, we manufacture components, medical devices, single-use devices, and tablets with active pharmaceutical ingredients, and we apply function-enhancing coatings to a number of products.

SP Medical is a company in SP Group A/S, an international plastics group headquartered in Denmark. SP Group is listed on NASDAQ OMX Copenhagen.

www.sp-medical.com